KOMBITEST™ FITC/PE Negative Control

Cat. No. / Kat. č.: ED7032

ENGLISH

1. Specification

Content	50 tests, 1 ml	
Usage	20 µl per test	
Specificity	Undefined epitope on a plant pathogen	
Clone	PPV-04	PPV-04
Isotype (mouse)	IgG2a	IgG2a
Fluorochrome	FITC	PE
λ excitation	488 nm	488 nm
Emission maximum	525 nm	575 nm

2. Intended use

The KOMBITEST™ FITC/PE Negative Control is designed for setting a boundary between negatively and positively stained lymphocytes when peripheral blood samples are analyzed using two-color (FITC/PE) KOMBITEST™ reagents.

3. Principle

Besides the specific binding of murine antibodies to the antigens, antibodies may also bind to human leukocytes through non-antigen (non-specific) means. This nonspecific binding of monoclonal antibodies may significantly aggravate the resolution between positive and negative fluorescent populations. By showing murine antibodies binding to human leukocytes through non-antigen (non-specific) means, the control sample enables establishing a boundary between negatively and positively stained leukocytes.

4. Specificity

The monoclonal antibody PPV-04 reacts with undefined epitope on a plant pathogen. The antigen is not present on human leukocytes. When added to human whole blood, the antibody PPV-04 does not bind specifically to antigens on the surface of human leukocytes.

5. Reagent provided

The reagent contains premixed combination of mouse monoclonal antibody PPV-04 labeled with Fluorescein isothiocyanate (FITC), and mouse monoclonal antibody PPV-04 labeled with R-phycoerythrin (PE). Labeled antibodies are diluted at optimum concentration in stabilizing phosphate

buffered saline (PBS) solution containing 15mM sodium azide. The content of a vial (1 ml) is sufficient for 50 tests.

6. Storage

Store vial in the dark at 2-8°C. Do not freeze.

7. Precautions

- Intended for In Vitro Diagnostic use in laboratories outside USA and Canada. This CE-IVD reagent is in conformity with the European In Vitro Diagnostic Medical Device Directive 98/79/EC.
- Stain the control sample under the same experimental conditions as the test sample. Use the same blood specimen.
- Do not use after expiration date stamped on vial label.
- Avoid prolonged exposure to light.
- The content of the vial must not freeze.
- Avoid contamination of the reagent.
- Any non-performance of staining protocol may produce false results.
- The reagent contains sodium azide (NaN₃) which is highly toxic in pure form. However, the concentration in the reagent (15mM) is not considered as hazardous. When disposing the reagent, flush the sink with large volume of water to avoid accumulation of explosive metal-azide in plumbing.
- Blood samples are considered as potentially infectious and must be handled with care. Avoid all contact of the sample with the skin, eyes and mucosa.

8. Necessary material not supplied

Material necessary for collection of peripheral blood, test tubes for staining of blood samples (e.g. 12 × 75 mm), automatic pipettes with disposable tips, vortex mixer, centrifuge, commercial lysing solution, phosphate buffered saline (PBS), flow cytometer.

9. Staining protocol

- Collect peripheral blood in a sterile tube with an anticoagulant (e.g. Heparin, EDTA).
- Add 20 µl of KOMBITEST™ FITC/PE Negative Control to a control tube.
- Add 100 µl of blood sample to the tube. Vortex the tube.

4. Incubate the tube for 20-30 minutes at room temperature in the dark.
5. Perform lysis of red cells using lysing solution. It is recommended to use a commercial lysing solution containing formaldehyde as a fixative. Follow the instructions of the lysing solution manufacturer.
6. Centrifuge the tube for 5 minutes at 300 g.
7. Remove supernatant and resuspend pellet with 3-4 ml of PBS.
8. Centrifuge the tube for 5 minutes at 300 g.
9. Remove supernatant and resuspend pellet with 0.3 – 0.5 ml of PBS.
10. Analyze the sample immediately using a flow cytometer or store sample at 2-8°C in the dark and analyze within 24 hours provided that cells were fixed.

10. Data analysis

Analyze the control sample stained with KOMBITEST™ FITC/PE Negative Control using a flow cytometer. Compensate fluorescent signals prior to or after data acquisition. Visualize compensated data in the side-scatter (SSC) versus forward-scatter (FSC) plot. Set the gate for lymphocyte population as shown in figure 1. Visualize lymphocytes in a dot-plot PE (FL2) versus FITC (FL1) as shown in figure 2. Set the gate corresponding to the negative population.

11. Limitations

- Flow cytometer may produce false results if the device has not been aligned and maintained appropriately.
- Data may be incorrectly interpreted if fluorescent signals were compensated wrongly or if gates were positioned inaccurately.
- Blood samples from abnormal patients may exhibit abnormal values of positive cells.
- Red blood cells from abnormal patients may be resistant to lysis using lysing solutions.
- In case of hyperleukocytose sample, it is recommended to dilute blood sample with PBS to obtain leukocyte density approximately 5×10^6 leukocytes/ml.
- Blood samples should be stained and analyzed within 24 hours from the blood collection.

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