



Reference: PTXBC000097

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BIGH3; CDB1; CDG2; CDGG1; CSD; CSD1; CSD2; CSD3; EBMD; LCD1; transforming growth factor-beta-induced protein ig-h3; RGD-CAP; RGD-containing collagen-associated protein; beta ig-h3; betaig-h3; kerato-epithelin; transforming growth factor beta-induced 68kDa		

For research use only.

GENE DESCRIPTION	transforming growth factor, beta-induced, 68kDa
GENE ID	7045
ORIGIN SPECIES	Human
ACCESSIONS	BC000097
VECTOR	pDONR223
NCBI SYMBOL	TGFBI
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr