

H2AFV-H2A histone family, member V Gene

Reference: PTXBC000098

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	H2A.Z-2; H2AV; histone H2A.V; H2A.F/Z; histone H2A.F/Z; purine-rich binding element protein B; H2A histone family member V		

For research use only.

GENE DESCRIPTION	H2A histone family, member V
GENE ID	94239
ORIGIN SPECIES	Human
ACCESSIONS	BC000098
VECTOR	pDONR223
NCBI SYMBOL	H2AFV
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctggaggcaagctggaaaggacagtgggaaggccaaggctaaaggcagatctcgtcacagagagctgggctacagtttctgtggccgcatccacagacactgaagactcgcaccacaagccatggaagggtgggtgccactgtgccgtgtacagtgtgcgattctgg
agttacctcactgcagagggtgctggagctggcaggtaatgcttaaggatctcaaaagtaaaagcgtatcactccgctcacttgcagcttgcattccgtgggtgatgaagagttggattctcttaaggtaccatagctgggggtggtgatccctcacatccacaatctctgattggaag
aaggacagcagaaaactgcttag