



Reference: PTXBC000119

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	proline-rich protein 14; proline rich 14		
For research use only.			

GENE DESCRIPTION	proline rich 14
GENE ID	78994
ORIGIN SPECIES	Human
ACCESSIONS	BC000119
VECTOR	pDONR223
NCBI SYMBOL	PRR14
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]