



Reference: PTXBC000125

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CED; DPD1; LAP; TGFB; TGFBeta; transforming growth factor beta-1; TGF-beta-1; latency-associated peptide; prepro-transforming growth factor beta-1; transforming growth factor beta 1		

For research use only.

GENE DESCRIPTION	transforming growth factor, beta 1
GENE ID	7040
ORIGIN SPECIES	Human
ACCESSIONS	BC000125
VECTOR	pDONR223
NCBI SYMBOL	TGFB1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr