

XPO5-exportin 5 Gene

Reference: PTXBC000129

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	exp5; exportin-5; ran-binding protein 21; exportin 5		
For research use only.			

GENE DESCRIPTION	exportin 5
GENE ID	57510
ORIGIN SPECIES	Human
ACCESSIONS	BC000129
VECTOR	pDONR223
NCBI SYMBOL	XPO5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaggctttcagaatggcaagttatcaacaaaggagcctgctgtgtggagaagatgaggctgcagatgaaacccagagctcaagagatgctggaggagcaactggtgaggatgttaaccgagaagtcagacctaatacaggtttgctgtttcaagaagggtgctgacc
acagtagtgcctccccagcagatggagacgatgaagaatgatggccacagaggtcacccctcagctatggcagagcttacagacctgggcaaatgtctgatgaagcatgaggatgtttacagcgtctattatcagcctcaatccctggcctggaagatactctgtcctccag
aggacaacctcacagctctgctgctcctcctcaaaagtgtgtcaggacactgctcgcagatgcagttacgtgtgtttcaccagtgctgctgaaaggcttacagatgcacgggcagcagcaggggtgcatggctccctggctccatctggcctccagatatacaggcactgcgccc
caggtacctggagataagagctgtatggagcaaatccctgaaalacagaaggactcactggaccagtttgactgcaagcttttaaacccctccctgcagaaagtggctgacaagcgcgaaaggaccaattcaaacgcctcattgctgttgctgattgggaacccctgggagagcagttc
cgaaaaagaagttcacattaagaatctccctcacttttcaaaaaaacaagccaatgctggagacggaggtgctggacaatgatgggggtgcctggccaccatctttgaacctga