



Reference: PTXBC000138

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CEAR; HNRNPM4; HNRPM; HNRPM4; HTGR1; NAGR1; hnRNP M; heterogeneous nuclear ribonucleoprotein M; CEA receptor; N-acetylglucosamine receptor 1; heterogenous nuclear ribonucleoprotein M4; hnRNA-binding protein M4		

For research use only.

GENE DESCRIPTION	heterogeneous nuclear ribonucleoprotein M
GENE ID	4670
ORIGIN SPECIES	Human
ACCESSIONS	BC000138
VECTOR	pDONR223
NCBI SYMBOL	HNRNPM
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr