

H1F0-H1 histone family, member 0 Gene

Reference: PTXBC000145

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	H10; H1FV; histone H1.0; H1.0, H1(0), H1-0; histone H1'; histone H1(0); H1 histone family member 0		

For research use only.

GENE DESCRIPTION	H1 histone family, member 0
GENE ID	3005
ORIGIN SPECIES	Human
ACCESSIONS	BC000145
VECTOR	pDONR223
NCBI SYMBOL	H1F0
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaccgagaattccagtcgccccctgcggccaagcccaagcgggccaagcctcgaagtcacagaccaccaagttatcagacatgatcgtggctccatccaggccgagaagaaccgcgtggtcctcgcgccagtcattcagaagtatatcaagagccactacaag
gtgggtgagaacgctgactcgagatcaagttgtccatcaagcgcctgtgaccaccggtgtcctcaagcagaccgaaggggtggggcctcgggtcctccgctagccaagagcgacgaaccaagaagtcagtggtcctcaagaagaccgaaggaatcaagaaggtagc
cacgccaagaaggtcatccaagccaagaaggtgctcctcaagcccaaccaagaacccaagccacccgggtcaagaaggccaagaagaagctgggtgccacgccaagaagcccaaaacccaagactgtcaagccaagccggtcaaggcatccaagcccaaaag
gccaaaccagtgaaccaccaagcaagtcagtgccaagggcgccaagaagaagtgga