

C1orf43-chromosome 1 open reading frame 43 Gene

Reference: PTXBC000152

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C1orf43; HSPC012; NICE-3; NS5ATP4; S863-3; HCV NS5A-transactivated protein 4; hepatitis C virus NS5A-transactivated protein 4; chromosome 1 open reading frame 43		

For research use only.

GENE DESCRIPTION	chromosome 1 open reading frame 43
GENE ID	25912
ORIGIN SPECIES	Human
ACCESSIONS	BC000152
VECTOR	pDONR223
NCBI SYMBOL	C1orf43
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgtccgcagtaactggctcctcggggggaatgtcgtgctggtgatggcctacgggagcctggtgttgactgctatttattttgtgaagaggcaaatatcgctttgcaatgaaatctcgaaggggacatgctcctgtgggacacaatgcccccaaggactgaaagaggaga
ttgatattcgactcaccagggttcaggatcaagatagagcccgctccttcagatgatgatgctagactactacaactggaaaccagggaatcaaatgttgcataactatctgtataggatgaaagctctggatgccattcgtaccttgagatccatttattctgaagccggcatcc
ccgttccttaattgggcaagaattccgctcctacgtcgtgatctgcgaaacactagtagcctttcaagggtgtacgcaaaagcactcattgatacccttttgatggctatgaacagccgctatgggacagggtctttggcagaatgagtacgtacgtacagaggccctgagtgag
ctggccactcgggttaagcacgaattgggagctcagcgacatcaccagtcagcagcaaaagacctaactcagtcctccctgaggtctcccaacaaccatccaggtgacatacctcccctcagtcagaagagtaaacgtgccaaagcacttcttgaaagagctttaaggataacta
taacacattggagagtactctgtga