

YWHAE-tyrosine 3-monooxygenase/tryptophan 5-monooxygenase activation protein, epsilon polypeptide Gene

Reference: PTXBC000179

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	14-3-3E; HEL2; KCIP-1; MDCR; MDS; 14-3-3 protein epsilon; 14-3-3 epsilon; epididymis luminal protein 2; mitochondrial import stimulation factor L subunit; protein kinase C inhibitor protein-1; tyrosine 3-monooxygenase/tryptophan 5-monooxygenase activation		

For research use only.

GENE DESCRIPTION	tyrosine 3-monooxygenase/tryptophan 5-monooxygenase activation protein, epsilon polypeptide
GENE ID	7531
ORIGIN SPECIES	Human
ACCESSIONS	BC000179
VECTOR	pDONR223
NCBI SYMBOL	YWHAE
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggatgatcgagaggatctgtgtaccagggcgaagctggccgagcagcgtgagcgatacgcgaatggtggagtcgaatgaagaagtagcaggggagtggtgagctgacagttgaagaagaaacccctctatctgtgcatataagaatgtgattggagctagaagagcctctg
gagaataatcagcagcattgaacagaagaagaacacagggaggagaagacaagctaaaaatgattcggaataatcgcaaatggtgagactgagctaaagttaattctgtgtgacattctggatgtactggacaaacacctcattccagcagctaacactggcaggtccaaggtttct
attataaaatgaaggaggactaccacaggtatctggcagaattgccacaggaacgacaggaaggagcgtgaggagaacagcctagtggtctataaagctgtagtgatattgcaatgacagaacttccaccaacgcacatctattcgttaggtcttctcattttccgtattctactacg
aaattcttaattccctgaccgtgctgcaggttgcaaaagcagctttgatgatgaattgcagaactggatacgcgtgagtgagaagactataaggactctacattatcatgcagttgttacgtgataatctgacactatgactcagacatgcagggtgacggtgaagagcagaataaa
gaagcgtgcagagcgtggaagacgaatcagtgga