

TMEM51-transmembrane protein 51 Gene

Reference: PTXBC000202

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C1orf72; transmembrane protein 51		
For research use only.			

GENE DESCRIPTION	transmembrane protein 51
GENE ID	55092
ORIGIN SPECIES	Human
ACCESSIONS	BC000202
VECTOR	pDONR223
NCBI SYMBOL	TMEM51
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatggccagtcacagccaatggctgcactatgcctgacgcctcggcctgggatgctgctcctggggtgatcatggccatgtggaacctgtaccggccttcagcgcggcgagagccaacagctcagggcagcaacaagaccgaggtgggtggcgcatcctcaag
agcaagaccttctctgtgacctacgtgctggctggggcgggggtgatgctgctcttctatctgctgagtaicagggataagaggaaagcagcggcaggcgagggacctggcccatgtccagcaccggacagggcctggcctcagccaggagggaagacagccaggagga
agaagaggaggatgaggagcctcctcaaggtactatgttccagctacgaggaagtgtgaacacaaactactcagaagcaaggagagagcagaaccggaggttgagcatctctccgtctctatgagtcactgacgggctcagcagaccacccacatccaccagg
ctgacgtggaggccagccctgggaacccctgacaggcagaacttaagtggccaaacgactgaaaccgctgaaagttcgaaggattaaatctgaaaagcttcacctcaagacttttaggatcaacctccagacaaaacgtccctctcctcagatagagcctttgactctccacc
gcagatgatgaagtcaggagaaggcccccacacccggccggccgactga