



Reference: PTXBC000212

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TFIIIC-BETA; TFIIIC110; general transcription factor 3C polypeptide 2; TF3C-beta; TFIIIC 110 kDa subunit; general transcription factor IIIC, polypeptide 2, beta 110kDa; transcription factor IIIC 110 kDa subunit; transcription factor IIIC subunit beta; gen		

For research use only.

GENE DESCRIPTION	general transcription factor IIIC, polypeptide 2, beta 110kDa
GENE ID	2976
ORIGIN SPECIES	Human
ACCESSIONS	BC000212
VECTOR	pDONR223
NCBI SYMBOL	GTF3C2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1