

RNF24-ring finger protein 24 Gene

Reference: PTXBC000213

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	G1L; RING finger protein 24		
For research use only.			

GENE DESCRIPTION	ring finger protein 24
GENE ID	11237
ORIGIN SPECIES	Human
ACCESSIONS	BC000213
VECTOR	pDONR223
NCBI SYMBOL	RNF24
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagctcggattccacattacaactcaggatgcctaattggattccagaatcgcctcacaatatattgtgggttttggtactgctatattgtcttcactccttagtttactctctgtgtgctacttgattaggctaagacatcaagcacacaagaattttatgcctacaacaggttatattaa
aagagaaagttaaagaattgaatttcatgagctctgtgcagtgctgctagaagactcaagcctcgagatgagttggggattggccatgtaagcacgcctccacagaaagtgccttattaagtggctggaggttcgtaaagtgtgtccctgtgcaacatgccagttctacagctggccca
gttcacagtaagcaggaccgtggacccctcagggggcccttcctggggcagagaacattgtatag