

C19orf43-chromosome 19 open reading frame 43 Gene

Reference: PTXBC000216

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C19orf43; fSAP18; My029 protein; chromosome 19 open reading frame 43		

For research use only.

GENE DESCRIPTION	chromosome 19 open reading frame 43
GENE ID	79002
ORIGIN SPECIES	Human
ACCESSIONS	BC000216
VECTOR	pDONR223
NCBI SYMBOL	C19orf43
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctccccgaggagacggcgaggcctcaggccgggaggtccgggccccggcggtggcggtggcgaggcgttgggctgagtcgggatcggggacgtcggcgagagcggggagggaggtgtggcgcggttcgagcccggttcggcgcggtga
actgttcgccaacgacggcagcttcctggagctgttcaagcggaagatggaggaggagcagcggcagcgagggagggagccggcggtccgcagcgacccgaccagtcggccgcccggctggccccgggatccgaagaggaggcggtccgggtccacacttag
cttcgtgggcaaacgcagagcggggaacaactagccctcaagacgggaatagtagccaagaagcagaagacgagggatgaggtattaacaagtaaaagtgacgcgtgggccaagtacatggcagaagtgaagaaagtaaaagctcaccagtcgggtgacgatataaaactcgg
ccccctggtgaaatga