

MIS12-MIS12, MIND kinetochore complex component, homolog (yeast) Gene

Reference: PTXBC000229

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MIS12, kinetochore complex component; homolog of yeast Mis12; MIS12, MIND kinetochore complex component, homolog; MIS12 homolog; protein MIS12 homolog; 2510025F08Rik; KNTC2AP; MTW1; hMis12		

For research use only.

GENE DESCRIPTION	MIS12, MIND kinetochore complex component, homolog (yeast)
GENE ID	79003
ORIGIN SPECIES	Human
ACCESSIONS	BC000229
VECTOR	pDONR223
NCBI SYMBOL	MIS12
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctgtgatccaatgacctacgagggccagttcttggcttcacgccacaaacgtgcattcggatctacattgcattcaagactacatttgaagtgtgagccgttgaaacaggttattctgaagaagctggatggcatccagactgtgacattagccagtcagattcgcaaa
tgcacagagaagttcttcttcattgaaaggacattttgataaccttttagcaaaatggagcaactgttttgcagctgattttacgtattccctcaaacatctgcttcctgaagataaatgtaaggagacaccttatagtgaggaaattttcagcatctccagaagaattgaacagttacagga
gaagtacaagactgaattatgtactaagcaggcccttctgcagaattagaagagcaaaaaattgtcaggccaaactcaaacagacgttgactttcttgatgagcttcataatgttggcagagatcatgggactagtgatttttagggagagtttagtatccctgttcagaactccagaaaacta
cagaacattagagacaatgtgaaaaggaaatcgaacgactgaaaatatcttaa