

DLK2-delta-like 2 homolog (Drosophila) Gene

Reference: PTXBC000230

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DLK-2; EGFL9; protein delta homolog 2; EGF-like protein 9; EGF-like-domain, multiple 9; delta-like 2 homolog; epidermal growth factor-like protein 9; delta like non-canonical Notch ligand 2		

For research use only.

GENE DESCRIPTION	delta-like 2 homolog (Drosophila)
GENE ID	65989
ORIGIN SPECIES	Human
ACCESSIONS	BC000230
VECTOR	pDONR223
NCBI SYMBOL	DLK2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgcttggctaacggtgccactgcctgacggcataaacgcttctcctgctctgagggcttggagcgttctgcacatcaacctggatgactgtgccagccgccatgccagagagggggcccgctgtcgggacgtgtccacgacttcgactgccttgcctcag
tggctatgggtggcaagacctgtgagcttcttacctgtccagaccccccaaccacagtggaacacctctagggccacacagctgtagtggtacctgccacggggccagccccccacagcgcaggggctggtctgtcgatctcagtgaaaggaggtgtcggaagcaagag
gctgggctaggtgagcctagcttggggccctggtgttggggccctcactgctgcccctggttctggctactgtgttgcctgacctgagggcctggcggcggtgtgtgtccccctggacctgtgtctacctgccccacactatgctccagctgccaggaccaggaggtgtcaggt
tagcatgtgccagcagggtccccctgccacgtgacttggccctgagcctggaagaccacagcactgtga