



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

MRPS7-mitochondrial ribosomal protein S7 Gene

Reference: PTXBC000241

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MRP-S; MRP-S7; RP-S7; RPMS7; S7mt; bMRP27a; 28S ribosomal protein S7, mitochondrial; 30S ribosomal protein S7 homolog; bMRP-27a; mitochondrial ribosomal protein S7		

For research use only.

GENE DESCRIPTION	mitochondrial ribosomal protein S7
GENE ID	51081
ORIGIN SPECIES	Human
ACCESSIONS	BC000241
VECTOR	pDONR223
NCBI SYMBOL	MRPS7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgtgtccccctcagatgaaggtgtcccgaggatgtcggcgctgcgttggcgctgcggcgggctgtcttcgacgttcacaggctaaactcagctgagatgagccgtatagtctgaattcaaggatcccttgattgacaaggaattatcgcgaagccagtgtaggagctaaactgag
gaggagaaatatgttcgggagctcaagaagactcagctcatcaaatgtctccgacagggaaaaacaaagtctgtgtttgaagcccgactcatcagtaaaatcccaacatgatgatgatagaggagaaacaaagtactgtgccagatccctcatgattcagactctggaagctgtgaaaaagga
agcagttgtgaagaattacatccgcttctcagagaagacaggcaaacatcgaaacgaacccctacacalcattccatcaagcattgaaaactctgtgacgtctatgtggcgctgtaccacatcctcaagggaaggccgtttcaccaggtccctgtaccctaccgacccggcgctccgcttc
ctagccatgaagtggtgatcatcctgagtgccgggataaaaagcaccagcggacactgatccggcgactgtcagacagctgtcggagctttccataacaggggcccgctgatcaagaggaagcatgacttgcaagaatgtgcagaggccaacctgtccctggcccataccgc
tggtgctgag

More info: Please contact us for any question: contact@proteogenix.fr