

MRPS23-mitochondrial ribosomal protein S23 Gene

Reference: PTXBC000242

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CGI-138; HSPC329; MRP-S23; 28S ribosomal protein S23, mitochondrial; S23mt; mitochondrial ribosomal protein S23		

For research use only.

GENE DESCRIPTION	mitochondrial ribosomal protein S23
GENE ID	51649
ORIGIN SPECIES	Human
ACCESSIONS	BC000242
VECTOR	pDONR223
NCBI SYMBOL	MRPS23
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcaggcagccgctggaaccgtaggagcatctctcggactcgggacgtgttcggcgccgggtgctgaaggagagccctgtggtttgacgtatatgacgctttcccccgtgaggagcccgcttccaaaggcctcgagtcgatatggcaagccaaagctccat
ccaagacatctggtaccagaggtatcggattagagcgaagttttatcaggtgtatgggtctgtcgaagagcctttgatctattcaatcacaactcaagctacgttcaacggtttgtggagaagtacactgagctacagaacttggagaacagatgaagagaagtatttggaaacagg
gaagcctttattggcagaaggtgtcattttaagacgagtaggcgaagcaaggactcaacacggaggtgacgtttcccggaatccgaacacttgaggtgtcagaccacagactgctgttggagaacacgagactcagaagaagtccacaggaccagcatttggaggcacctgcag
accagtcgaaggtctcttgcctccctga