

NME3-non-metastatic cells 3, protein expressed in Gene

Reference: PTXBC000250

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DR-nm23; NDPK-C; NDPKC; NM23-H3; NM23H3; c371H6.2; nucleoside diphosphate kinase 3; NDK 3; NDP kinase 3; NDP kinase C; non-metastatic cells 3, protein expressed in; nucleoside diphosphate kinase C; NME/NM23 nucleoside diphosphate kinase 3		

For research use only.

GENE DESCRIPTION	non-metastatic cells 3, protein expressed in
GENE ID	4832
ORIGIN SPECIES	Human
ACCESSIONS	BC000250
VECTOR	pDONR223
NCBI SYMBOL	NME3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatcgtcctggtgctgaccatcttcgtaacctctccccgcggcctgcaccggcgcaacgaacgcaccttcctggcctgaagccgggacggcgtgcagcggcggtgggtggcgagattgtgcggcgcttcgagaggaagggttcaagttgggtggcgtgaagctggtgcagg
cctccgaggagctgctgctgagcactacgccgagctgcgtgaacggcgttctacggccgctgtcaagtatatggcctccggcggtgggtggcctggtttggcaggggctggacgtgggtgcgacctcgcggcgctcatcgagccagcaacccggccgacgccccgccc
ggcaccatccgggggatttctgcatcgaggttggaagaacctgattcacggcagcgactcgggtggagagtgcccgcgcgagatcgtctctgttccgcgcagacgagctcctctgctgggaggaacagcgtgggcactggctgtatgagtag