

NDUFA1-NADH dehydrogenase (ubiquinone) 1 alpha subcomplex, 1, 7.5kDa Gene

Reference: PTXBC000266

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CI-MWFE; MWFE; ZNF183; NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 1; NADH dehydrogenase (ubiquinone) 1 alpha subcomplex, 1, 7.5kDa; NADH oxidoreductase subunit MWFE; NADH-ubiquinone oxidoreductase MWFE subunit; NADH:ubiquinone oxidoreductase		

For research use only.

GENE DESCRIPTION	NADH dehydrogenase (ubiquinone) 1 alpha subcomplex, 1, 7.5kDa
GENE ID	4694
ORIGIN SPECIES	Human
ACCESSIONS	BC000266
VECTOR	pDONR223
NCBI SYMBOL	NDUFA1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtggttcgagattccccggactctcgtcatggcggtgtgtgtgtgattccaggactgctactgcgtacatccacaggttcactaacggggcaaggaaaaagggtgtctcatttgggtatcactggagctctgatgaaagagataggcgcatctctggagttgatcgttactatgtgtcaaaagggttggagaacattgattaa