

GTF2A2-general transcription factor IIA, 2, 12kDa Gene

Reference: PTXBC000287

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HsT18745; T18745; TF2A2; TFIIA; TFIIA-12; TFIIA-gamma; TFIIAS; transcription initiation factor IIA subunit 2; TFIIA gamma subunit; TFIIA p12 subunit; general transcription factor IIA 2; general transcription factor IIA, 2, 12kDa; transcription initiation		

For research use only.

GENE DESCRIPTION	general transcription factor IIA, 2, 12kDa
GENE ID	2958
ORIGIN SPECIES	Human
ACCESSIONS	BC000287
VECTOR	pDONR223
NCBI SYMBOL	GTF2A2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcatatcagttatcacagaataactacttgggaaacagctctcaggagagcctagatgagctcatcagctctcaacagatcaccccccaacttgccctcaagttctactcagtttgataaggctataaatgcagcactggctcagagggtcagggaacagagtcatttcagggtctctctaaatcgtacagattctgcgataatgtgtggactttgtactgaatgatgtgaattcagagaggtgacagaacttataaagtgagataaagtgaataattgtagcctgtgatggtataaatactggctccaatactacagaatga