

### KLHDC3-kelch domain containing 3 Gene

Reference: PTXBC000295

|                |                                                                                                                     |                     |                |
|----------------|---------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                 | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | PEAS; dJ20C7.3; kelch domain-containing protein 3; testis intracellular mediator protein; kelch domain containing 3 |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | kelch domain containing 3                                                                         |
| <b>GENE ID</b>          | 116138                                                                                            |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC000295                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | KLHDC3                                                                                            |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

#### ORF (Only for reference, product provided without STOP codon )

atgttacggtgacagctgcacgtggaggcgggccccgcagggtgaaccatgctgcagtgctgtcgggcatcgggtatactccttgggggttactgctctgtgaagactatgagacactgcgtcagatagatgtgcacatttcaatgcagtgctcttgcgttgacaaagctgcccc  
ggfgaagtctgccatcctgtggcaagctcctgtgtaccctacatgcctatggacactcaaccgtctctcagcagacacagctcctcttggggcggggaatgacaccgaaggccctgcaatgtgctctatgaccttgcataacgcacaagtggttcacacccgagtgctcag  
ggacagttctggggccccgggatggacattcagcctgtgtcctaggcaagatcatgtacatttttgggggctacgagcagcaggcggactgttttccaatgacattcacaagctagataccagcaccaatgacatggactcttctgtacaaaggcagccctgcacgctggaggacttc  
actcagcccaatgctgggaagtcacatgtatgtcttgggggcccgtgcgcgaccgtttgggcccattcccaacaatgagatttactgcaaccgattcagcttggacaccagaactgaggcttggctggactgtccccgactccagtgctgctgaggggcccggagccactgc  
gcttggctacaatggggagctgtacatcttgggtgtataatgaaggctgaaccggcacttccatgacctctggaagttaactcctgtgtcttacctgaaaaagattgaaccgaagggaaggccatgtccccgcgcgcccagtgctgtctgtattgttgacaagattgtcct  
cttgggggtaccagtcctcctgaggaaggcctgggagatgaattgaccttatagatcattctgacttacacatttggactttagccctagctctgaagactctgtgcaaaactggccgtgattcagtagataacctagaccagtcctgttgcctcatgatatcagctgggagctgaatgccatga  
ccaccaacagcaatatcagtcgccccatcgtctcctccatgggtag