



Reference: PTXBC000296

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DFNA67; DNFA67; ORP-2; ORP2; oxysterol-binding protein-related protein 2; OSBP-related protein 2; oxysterol binding protein like 2		

For research use only.

GENE DESCRIPTION	oxysterol binding protein-like 2
GENE ID	9885
ORIGIN SPECIES	Human
ACCESSIONS	BC000296
VECTOR	pDONR223
NCBI SYMBOL	OSBPL2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1