

RTN1-reticulon 1 Gene

Reference: PTXBC000314

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NSP; reticulon-1; neuroendocrine-specific protein; reticulon 1		
For research use only.			

GENE DESCRIPTION	reticulon 1
GENE ID	6252
ORIGIN SPECIES	Human
ACCESSIONS	BC000314
VECTOR	pDONR223
NCBI SYMBOL	RTN1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggactgtgtgtggagcaactggaaaagtcaggctattgacctgtgtattggcgggacatcaagcagacgggcatcgtgtttgggagtttctcgtctgtctcttctccctgacccagttcagcgtgtgtgagcgtcgtggcctacctggccctggccgcactctcagccaccatcagtttcc
gcatctacaagtcgtgtttacaagcagtcgagaaaaccgacgaagggccacccttcaaggcctacttggagcttgagatcacccttctcaggagcagattcagaagtacacggactgcctgcagttctacgtgaacagcacacttaagggaactgaggagctcttcttgcaggacctgg
tggattccttaaaattgagtcctgatgtggctcctgacctacgttggcgtctcttcaatggcctgacctgctctcatggctgtgtgttcaatgtttactctacctgtatgtatgttaagcaccaggcacagattgaccaatatctgggactgtgaggaactcacataaatgtgtgtggcaaa
gattcaggctaaaatccaggcgctaagaggcacgctgagtaa