

DCTN3-dynactin 3 (p22) Gene

Reference: PTXBC000319

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DCTN-22; DCTN22; dynactin subunit 3; dynactin 3 (p22); dynactin complex subunit 22 kDa subunit; dynactin light chain		

For research use only.

GENE DESCRIPTION	dynactin 3 (p22)
GENE ID	11258
ORIGIN SPECIES	Human
ACCESSIONS	BC000319
VECTOR	pDONR223
NCBI SYMBOL	DCTN3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgggtctgactgacttcagcgggtacagcccgagtggaagagctggagcgtgggtgtacggccggggggcgcgcggtcacggaaggtggctgacggcctgggtcaaggtgcaggtggcttggggaacattccagcaagaggagaggggtgaagattcttaca
aaaagattgaagatctgatcaagtaactggatcctgagtacatcgaaccgattgccatactgatgccttaagctgcaattcatcctagcagaggagcagtttatccttcccaggttgactcctggagcaggtgaatgccttgggtcccatgctggacagtgtcacatcaaacccgttct
gagcatgctgcccgctgacggttggccagatccacattcagcagcaggtctcatgggagtgaggagtcctgatgaggcaggaagtttagtgaagatgtggccttggccagttccttctgtgctacacttggccctacaggaccagtgtgtggaatcactga