

SFN-stratifin Gene

Reference: PTXBC000329

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	YWHAS; 14-3-3 protein sigma; 14-3-3 sigma; epithelial cell marker protein 1		
For research use only.			

GENE DESCRIPTION	stratifin
GENE ID	2810
ORIGIN SPECIES	Human
ACCESSIONS	BC000329
VECTOR	pDONR223
NCBI SYMBOL	SFN
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagagagccagctctgatccagaaggccaagctggcagagcagggccgaacgctatgaggacatggcagcctcatgaaaggcgcgtggagaaggccgaggagctctctgcgaagagcgaacctgctctcagtagcctataagaacgtggfggggccagagggtgcct
ggagggtgctgtccagtattgagcagaaaaaacaacgaggagggctcggaggagaaggggcccgagggtgcgtgagtagccggagaagggtggagactgagctccagggtgctgacacccgtgctgggctgctggacagccacatcaaggaggccgggacgccgagag
ccgggtctctacctgaagatgaagggtgactactacgtacgtggccgaggtggccaccgggtgacgacaagaagcgcatcattgactcagccgggtcagcctaccaggaggccatggacatcagcaagaaggagatgcccccacaccccatccgctgggctggcctgaa
cttttcgcttctcactacgagatcgccaacagccccgaggggccatctctgccaagaccatttcgacgaggccatggctgatctgcacacctcagcgaggactcctacaagacagaccctcatcagctgctgagagacaacctgacactgtggacggccgacaacgc
cggggaagagggggcgagggtccccaggagccccagagctga