

## RPLP0-ribosomal protein, large, P0 Gene

Reference: PTXBC000345

|                |                                                                                                                                                                                                                                        |                     |                |
|----------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                                                    | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                                                   | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | L10E; LP0; PRLP0; RPP0; 60S acidic ribosomal protein P0; 60S ribosomal protein L10E; acidic ribosomal phosphoprotein P0; neutral ribosomal phosphoprotein P0; ribosomal protein, large, P0; ribosomal protein lateral stalk subunit P0 |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | ribosomal protein, large, P0                                                                      |
| <b>GENE ID</b>          | 6175                                                                                              |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC000345                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | RPLP0                                                                                             |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atgccagggaagacagggcgacctggaagtccaactcttcaagatcatccaactattggtgattatccgaaatgtttcattgtgggagcagacaatgtgggtccaagcagatgcagcagatccgcatgtccctcgtgggaaggctgtggtgctgatgggcaagaacaccatgat  
gcgcaaggccatccgagggcacttggaatacaaccagctctggagaactctgctctatccgggggaatgtgggtcttctgttcccaaggagacactcactgagatcagggacatgttctggtcccaataaggtgccagctgctgcccgtgtggtgccattgccccatgtgaagt  
cactgtgccagcccagaacactggtctcgggccgagaagacctcttttccaggctttaggtatcaccactaaaatctccaggggaccattgaaatcctgagtgatgtgcagctgatcaagactggagacaaagtgggagccagcgaagccagctgctgaacatgctcaacatctcc  
cccttctcttgggtggtcgcagcagctgttcgacaatggcagcatctacaacctgaagtgttgatcacagaggaaactctgcatctcgtctcctggagggtgtccgcaatgttgcagtgctgtctgcagattggctacccaactgttgcacgtacccattctatcatcaacg  
gggtacaacgagtcctggtctgtgtggagacgattacaccttccactctgctgaaaaggtcaaggccttcttggctgatccatctgcttctgtgctgctgcccctgtggtgctgctgccaccacagctgctcctgctgctgctgcagcccagctaaaggtgaagccaagggaagtcgg  
aggagtcgacgaggatattggatttggctcttgaactaa