

EIF4E2-eukaryotic translation initiation factor 4E family member 2 Gene

Reference: PTXBC000360

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	4E-LP; 4EHP; EIF4EL3; IF4e; eukaryotic translation initiation factor 4E type 2; eIF-4E type 2; eIF4E-like cap-binding protein; eIF4E-like protein 4E-LP; eukaryotic translation initiation factor 4E homologous protein; eukaryotic translation initiation fact		

For research use only.

GENE DESCRIPTION	eukaryotic translation initiation factor 4E family member 2
GENE ID	9470
ORIGIN SPECIES	Human
ACCESSIONS	BC000360
VECTOR	pDONR223
NCBI SYMBOL	EIF4E2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaacaacaagttcgacgtttgaaagatgatgacagtgaggaccatgacagaatgaagaaaacagcacagaaagatggtgagaaggaaaaacggaacgagacaagaatcagagcagtagcaagagaaggctgtgtccctggaccggcagagcatccctgcagtacaa
ctacacttttgggtactccaggagaacccccggccgtccacgagctcagagctatgaacagaatafcaaacagattggcactttgctctgtggagcagttctggaggttttatagccacatggtacgtctctggggacctgacaggccacagtgacttccatctctcaagaagggaatt
aaacccatgtgggaggtgatgcaataaaaatgggtgcaagtggtatttggctcggaaggccttggcctccgttgcctggagaatctcatgttggccatgctgggggaacatggttggggagagatctgtgggctgtggtctgtccgtttcaggagacattattcaa
tatggataagactgccagtgaacgaagcaaccacagccgaatccgggacacacttcggcgagtgcttaacctacctcccaaccattatggaatacaaaactcacaccgacagcatcaaatgccaggcagctgggccccaaaggctcttttcaaaacctctggaagccggg
ttgaatgtccatgacctctccctctctggatggcaccatcattga