

RPL18-ribosomal protein L18 Gene

Reference: PTXBC000374

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	L18; 60S ribosomal protein L18		

For research use only.

GENE DESCRIPTION	ribosomal protein L18
GENE ID	6141
ORIGIN SPECIES	Human
ACCESSIONS	BC000374
VECTOR	pDONR223
NCBI SYMBOL	RPL18
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggagtgacatccgccataacaaggaccgaaaggttcggcgcaaggagcccaagagccagatctacgtgaggtgttggtcaagttatacaggtttctggcagaagaaccaactccacattcaaccagggtgtgtgaagaggtgtttatgagtcgaccaaccggccgctc
tgtcccttcccgatgatccggaagatgaagcttctggccgggaaaacaagacggcgtgtgttggtgggaccataactgatgatgtgcgggttcaggaggtaccctgaaggtatgtgcactgcgcgtgaccagccggcccgagccgcacatctcaggcagggggcaagat
cctcacttgcaccagctggccctggactccccctaagggctgtggcactgtctgtctcgggtcctcgcaaggccgagaggtgtacggcatttcggcaaggcccccagggaacccgcacagccacacaaacctacgtcgcctcaaggccgggaagttcgagcgtgccagagg
ccgacggggcagccgaggtacaaaaactaa