

WTAP-Wilms tumor 1 associated protein Gene

Reference: PTXBC000383

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	pre-mRNA-splicing regulator WTAP; Mum2; PNAS-132; WT1-associated protein; Wilms' tumour 1-associating protein; female-lethal(2)D homolog; hFL(2)D; wilms tumor 1-associating protein; Wilms tumor 1 associated protein		

For research use only.

GENE DESCRIPTION	Wilms tumor 1 associated protein
GENE ID	9589
ORIGIN SPECIES	Human
ACCESSIONS	BC000383
VECTOR	pDONR223
NCBI SYMBOL	WTAP
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaccaacgaagaacctctccaagaagggtcgaatgagtgaaacagacttcaagttatggcaagagatgagtttaattctaatggaacaatatgaagcatatgtacaagcttggagggaagtaacacagatcttaacttaataatgatgtaactggcctaagagagtgtaagaaaaac
taagcaacaacagcaggagtgctgcagcagggaaaacatcttgtatgcgactagcaaccaaggaacaagagatgcaagagtgactactcaaatccagtaacctcaagcaagtcagcagccgagcgttgcccaactgagatcaacaatggttagcccgatcaactgtttttct
aaaatgaaagggaactggaacagactaagacaactggaacaagcccaaatgaactgagtgctggaagtttacgcctgataggtaa