

TSPAN4-tetraspanin 4 Gene

Reference: PTXBC000389

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NAG-2; NAG2; TETRASPAN; TM4SF7; TSPAN-4; tetraspanin-4; novel antigen 2; tetraspan TM4SF; transmembrane 4 superfamily member 7; tetraspanin 4		

For research use only.

GENE DESCRIPTION	tetraspanin 4
GENE ID	7106
ORIGIN SPECIES	Human
ACCESSIONS	BC000389
VECTOR	pDONR223
NCBI SYMBOL	TSPAN4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgcgcgctgctccaggccgtcaagtacatgttcgcttcaacctgctcttctggctggaggctgtggcgtgctgggtgtcggcatctgctggccgccacacagggagcttcgccagctgtccttcttccctgctcgtgctgccaactgctcatcaccggc
gccttctcatggccatcggtctgtggctgctgggtgccatcaaggagaacaagtgcctcctgctcactttcttctgctgctgctggtgttctgctgagggccaccatcgccatcctcttcttgcctacacggacaagattgacaggtatgccagcaagacctgaagaaggctt
gcacctgtacggcagcagggcaactgtggcctcaccacgcctggagcatcatcagaccgacttccgctgctgtggcgtctccaactacactgactggttcgaggtgtacaacgccacgagggtacctgactcctgctgtggagttcagtgagagctgtgggtgcacgccccg
gcacctgtggaaggcgccgtgtctacgagacgggtgaagggtgtggcttcaggagaacctgctggctgtgggcatcttgggctgtgcacggcgctggtgcagatcctgggctgaccttcgcatgacctgtactccaagtgtgtaaggcagacacctactgcgcgtag