

FAM127B-family with sequence similarity 127, member B Gene

Reference: PTXBC000393

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein FAM127B; CXX1b; MAR8A; mammalian retrotransposon derived protein 8A; family with sequence similarity 127 member B		

For research use only.

GENE DESCRIPTION	family with sequence similarity 127, member B
GENE ID	26071
ORIGIN SPECIES	Human
ACCESSIONS	BC000393
VECTOR	pDONR223
NCBI SYMBOL	FAM127B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacgggtcgggtgcagctgatgaaggccctctggccgggccccctccggcccgccgctcgtcgtgaggaacccgattcccttccccgagacgtttgacggagataccgaccgactcccgaggttcacgtgcagacgtgctcctacatgttcgtggacgagaacacgttctcaa
cgacgccctgaagtgacgttctcatcacccgctcaggggccagccctgcagtggtgatccctacatcaggaaggagagccccctgctcaatgattaccggggcttctcgtgccgagatgaagcgggtcttggatggaggaggacgaggacttctag