



Reference: PTXBC000397

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GPATC5; GPATCH5; Spp2; T54; G patch domain and KOW motifs-containing protein; protein MOS2 homolog; G-patch domain and KOW motifs		

For research use only.

GENE DESCRIPTION	G patch domain and KOW motifs
GENE ID	27238
ORIGIN SPECIES	Human
ACCESSIONS	BC000397
VECTOR	pDONR223
NCBI SYMBOL	GPKOW
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]