

PAFAH1B2-platelet-activating factor acetylhydrolase, isoform Ib, beta subunit 30kDa Gene

Reference: PTXBC000398

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HEL-S-303; platelet-activating factor acetylhydrolase IB subunit beta; PAF acetylhydrolase 30 kDa subunit; PAF-AH 30 kDa subunit; PAF-AH subunit beta; PAF-AH1b alpha 2 subunit; PAFAH subunit beta; epididymis secretory protein Li 303; intracellular platele		

For research use only.

GENE DESCRIPTION	platelet-activating factor acetylhydrolase, isoform Ib, beta subunit 30kDa
GENE ID	5049
ORIGIN SPECIES	Human
ACCESSIONS	BC000398
VECTOR	pDONR223
NCBI SYMBOL	PAFAH1B2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagccaaggagactcaaacccagcagctattccgcatgcagcagaagatattcaaggagatgacgatggatgtctcagcacaacagattgtttggactgtaaagacaagagcctgatgtactgttcgtgggagactccatggtgcagttaatgcagcaatagagatggcgaga
gctttttcccacttcactgaatttgaattggggagatacaacaagacatgtttgtggagactaaagaatggagaactggagaatattaagcctaaggctcattgttctgtggtaggaacaataaccacgaaaatacagcagaagaagtagcaggtgggatcgaggccattgta
caacttatcaacacaaggcagccacaggccaaatcattgtattgggtttgtacctcaggtgagaaacccaatcctttgaggcaaaagaacgccaaggtgaaccaactcctcaaggtttcgctgccgaagcttgccaacgtgcagctcctggataccgacgggggtttgtgcactcgga
cgggtccatctctgccacgacatgtttgattttctgcatctgacaggagggggctatgcaaaagatctgcaaacccctgcatgaactgatcatgcagttgttggaggaaacacctgaggagaacaaccacattgcctga