

## MAPK13-mitogen-activated protein kinase 13 Gene

Reference: PTXBC000433

|                |                                                                                                                                                                                                    |                     |                |
|----------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                               | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | MAPK 13; MAPK-13; PRKM13; SAPK4; p38delta; mitogen-activated protein kinase 13; MAP kinase 13; MAP kinase p38 delta; mitogen-activated protein kinase p38 delta; stress-activated protein kinase 4 |                     |                |

**For research use only.**

---

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | mitogen-activated protein kinase 13                                                              |
| <b>GENE ID</b>          | 5603                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC000433                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | MAPK13                                                                                           |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atgagcctcatcggaaaaaggcttctacaagcaggacgtcaacaagaccgctgggagctgcccagacacgtcggcagcggggcctatggctcgtgtgtcggccatcgacaagcggtcaggggagaagggtggccatcaagaagctgagccgac  
cccttcagtcaggatcttcgccaagcgcgcctacgggagctgctgctgctgaagcacatgcagcatgagaacgtcattgggctcctggtatgtcttccccagcctcctcctgcgcaacttctatgacttctactcgtgtgatgccctcatgcagacggatctgcagaagatcatgggat  
ggagttcagtgaggagaagatccagtacctggtgtatcagatgctcaaaaggccttaagtacatccactctgctgggctcgtgcacagggacctgaagccaggcaacctggtgtgaatgaggactgtgaactgaagattctggatttgggctggcgcacatgcagacccgagatgact  
ggctacgtggtgacccgctgggtaccgagcccccaggtgatcctcagctggatgcactacaaccagacagtgagacatctggtctgtgggctgtatcatggcagagatgctgacagggaaaactctgttcaaggggaaagattacctggaccagctgaccagatcctgaaagtgaccgg  
gggtcctggcagggagtttgtgcagaagctgaacgacaaagcggcctaatcctacatccagtcctgccagacccccagggaaggatttactcagctgttcccacgggcccagccccaggctgcggacctgctggagaagatgctggagctagacgtggacaagcctgacggc  
cgcgagggcctcaccatccctctttgaacccttcgggacccctgagggaagagcggagggccagcagccgtttgatgattcctagaacacgagaactcagatggatgaatggaagcagcacatctacaaggagattgtgaactcagccccattgccgggaaggactcagggc  
ggcgagtggtgcatgaagctgttag