

ENSA-endosulfine alpha Gene

Reference: PTXBC000436

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ARPP-19e; alpha-endosulfine; endosulfine-alpha variant 1; endosulfine-alpha variant 3; endosulfine alpha		

For research use only.

GENE DESCRIPTION	endosulfine alpha
GENE ID	2029
ORIGIN SPECIES	Human
ACCESSIONS	BC000436
VECTOR	pDONR223
NCBI SYMBOL	ENSA
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtccagaacaagaagaaccctgcggaggagaccggcgaggagaagcaggacacgcaggagagaagaagggtattctgcctgagagagctgaagaggcaagctaaaggccaatacccaagcctaggacaaaagcctggaggtccgacttctcatgaagagactc
cagaaaggcgaagactttgactcaggagactacaacatggcgaagccaagatgaagaataagcagctgccaaagtcaggaccagacaagaacctggtgactggtgatcacatccccacccacagatctgccccagagaaagtcctcgtctcaccagcaagcttcgggt
ggccaagtgaatga