

RPS13-ribosomal protein S13 Gene

Reference: PTXBC000475

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	S13; 40S ribosomal protein S13; ribosomal protein S13		
For research use only.			

GENE DESCRIPTION	ribosomal protein S13
GENE ID	6207
ORIGIN SPECIES	Human
ACCESSIONS	BC000475
VECTOR	pDONR223
NCBI SYMBOL	RPS13
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggtcgcacatgctccgggaaggcctgtccagtcgggtttaccctatcgacgcagcgtcccttgggtgaagtgacatcgacgacgtgaaggagcagattacaactggccaagaaggccttactcctcacagatcgggtgaatcctgagagattcacatgggttgca
caagtacgttttgacaggcaataaaatttaagaattcttaagctaaaggactgtctctgactctcctgaagatctctaccatttaataagaagcagttgctgtcgaagcatctgagagggaacagaaaggataaggatgctaaattccgtctgattctaataagagccggattcacctg
ttggctcgatattataagaccaagcagtcctccctccaattggaatatgaatcatctacagcctctgcccctggcgcataa