

RPL9-ribosomal protein L9 Gene

Reference: PTXBC000483

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NPC-A-16; 60S ribosomal protein L9		
For research use only.			

GENE DESCRIPTION	ribosomal protein L9
GENE ID	6133
ORIGIN SPECIES	Human
ACCESSIONS	BC000483
VECTOR	pDONR223
NCBI SYMBOL	RPL9
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagactattctcagcaatcagactgtcgacattccagaaatgtcgacattactctgaaggacgcacagttatcgtgaaggcccccagaggaaacctgcggagggaactcaatcacatcaatgtagaactcagcctcttggaagaaaaaaagaggctccgggtgacaaatgt
gggtaacagaaaggaaactgctaccgttcggactattttagtcatgtacagaacatgatcaagggtgttacactgggctccgttacaagatgaggtctgtgtatgctcaacttccccatcaacgttgttatccaggagaatgggtctctgtgaaatccgaaattcttgggtgaaaaatacatc
cgagggttcggatgagaccaggtgtgtgtgtcagttatctcaagcccagaaagatgaattaatccttgaaggaaatgacattgagctgtttcaaatcagcggccttgattcagcaagccacaacagttaaaaacaaggatcaggaatttttgatggtatctatgtctgaaaaaggaa
ctgttcagcaggctgatgaataa