

SNRPD2-small nuclear ribonucleoprotein D2 polypeptide 16.5kDa Gene

Reference: PTXBC000486

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SMD2; SNRPD1; Sm-D2; small nuclear ribonucleoprotein Sm D2; small nuclear ribonucleoprotein D2 polypeptide 16.5kDa; snRNP core protein D2; small nuclear ribonucleoprotein D2 polypeptide		

For research use only.

GENE DESCRIPTION	small nuclear ribonucleoprotein D2 polypeptide 16.5kDa
GENE ID	6633
ORIGIN SPECIES	Human
ACCESSIONS	BC000486
VECTOR	pDONR223
NCBI SYMBOL	SNRPD2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagcctctcaacaagccaagagtgagatgacccagaggagctgcagaagcgagaggagggaatttaacaccgggtccactctctgtgtcacacagtcagtcagaacaataccaagtgctcatcaactgccgcaacaataagaaactcctggcgcggtgaaggccttcg
ataggcactgcaacatgggtctggagaacgtgaaggagatgtggactgaggtacccaagatggcaaggccaagaagaagtcgaagccagtcacaaagaccgctacatctccaagatgttcctgcgcggggactcagtcctgtgtcctgcggaaaccgctcatcgccggcaagt
ag