



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

MGST3-microsomal glutathione S-transferase 3 Gene

Reference: PTXBC000505

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GST-III; microsomal glutathione S-transferase 3; microsomal GST-3; microsomal GST-III; microsomal glutathione S-transferase III		

For research use only.

GENE DESCRIPTION	microsomal glutathione S-transferase 3
GENE ID	4259
ORIGIN SPECIES	Human
ACCESSIONS	BC000505
VECTOR	pDONR223
NCBI SYMBOL	MGST3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgtccctcctaaggaatatgtgtttgtgcttctaactggctgctgcagcgttataatggtgccaccctagccatcaatgttccaaggcccgcaagaagtacaaagtgaggatctctatcatgtacagcagcgaccctgaaaaatgggcacatcttcaactgattcacgagcccccag
aacacgttgcaagaagtatcccctctcttattttctatgttgagggtgtttacacccctgtatagctctgtcgtcgctgctgcgtgattgttgagacagcttcttatgtctatgcacgggaagcaacggcaagcgtatgctcgaggaagccctgggtccatcgccctctgggctt
gtggggcacaacatctgtctctctctgcacacatctgtgttggctttaaaggtgtggcgatggagcccaaatctgtccataaa