



Reference: PTXBC000513

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LNCR2; NACHRA3; PAOD2; neuronal acetylcholine receptor subunit alpha-3; cholinergic receptor, nicotinic alpha 3; cholinergic receptor, nicotinic, alpha 3 (neuronal); cholinergic receptor, nicotinic, alpha polypeptide 3; neuronal nicotinic acetylcholine re		

GENE DESCRIPTION	cholinergic receptor, nicotinic, alpha 3
GENE ID	1136
ORIGIN SPECIES	Human
ACCESSIONS	BC000513
VECTOR	pDONR223
NCBI SYMBOL	CHRNA3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1