

RPS24-ribosomal protein S24 Gene

Reference: PTXBC000523

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DBA3; S24; 40S ribosomal protein S24		

For research use only.

GENE DESCRIPTION	ribosomal protein S24
GENE ID	6229
ORIGIN SPECIES	Human
ACCESSIONS	BC000523
VECTOR	pDONR223
NCBI SYMBOL	RPS24
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaacgacacgtaactatccgactagaaagttcatgaccaaccgactactcagaggaaacaatggtcattgatgtcttcacccgggaaggcgacagtcgctaagacagaaatcgggaaaaactagccaaatgtacaagaccacacggatgtcatcttgtatttgattcaga
actcattttgggtggcgaagacaactggcttggcatgatttatgattccctggattatgcaagaaaaatgaacccaacataacttgcaagacatggcctgtatgagaagaaaaagacctcaagaagcaacgaaaggaacgcaagaacagaatgaagaaagtcaggggactgca
aaggccaatgtgtgctgcaaaaagtga