

TMEM98-transmembrane protein 98 Gene

Reference: PTXBC000526

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TADA1; transmembrane protein 98		
For research use only.			

GENE DESCRIPTION	transmembrane protein 98
GENE ID	26022
ORIGIN SPECIES	Human
ACCESSIONS	BC000526
VECTOR	pDONR223
NCBI SYMBOL	TMEM98
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagactgtggtgattgtgcataggtgtgctggccaccatcttctggcttctgttgcagccttgggtgttgcaggcagcgtactgccggccgcgagacctgctgcagcgtatgattctaagccattgtggacctcattgtgcatggagaccagctgtgagccctctgagtttgaacatggacgatgtcttatcaccaacccccacattgaggccattctggagaatgaagactggatcgaagatgctcgggtctcatgtccactgcatgcatcttgaagattgtcacactgtgacagagaagctgttggcatgacaatgggtcttgggccaagatgaagactcagcca
gtgtcagcgacatcattgtgtggtgccaagcggatcagccccaggggtggtggtggaagtcgatgtaccctcgttggaccccaaacctcgtgacgcacgacgactgccctgctctgtgtgtcagtcacctgtgtgtgtgacaaggaatgcctgcatctgacgggagggcctgga
ctggattgaccagctctgtcgtgctgaggagcatttgaagtccttcgagaagcagccctagcttctgagccagataaaggcctccagggcctgaaggcttctgcaggagcagctctgcaatttag