

RPL19-ribosomal protein L19 Gene

Reference: PTXBC000530

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	L19; 60S ribosomal protein L19; ribosomal protein L19		
For research use only.			

GENE DESCRIPTION	ribosomal protein L19
GENE ID	6143
ORIGIN SPECIES	Human
ACCESSIONS	BC000530
VECTOR	pDONR223
NCBI SYMBOL	RPL19
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagtatgctcaggcttcagaagaggctcctctagtgctcctcgtgtggcaagaaggtctgttagacccaatgagaccaatgaatcgccaatgccactccgtcagcagatccggaagctcaaaagatgggtgatccgcaagcctgtgacgttcattccggg
ctcgaatgccggaataacaccttggcccgccggaaggcaggcacaatggcataggttaagcggaggttacagccaatgccgaatgccagagaaggtcacatggatgaggagaatgaggaatttgcgccgctgtcagaagataccgtgaatctaagaagatgatgccacatgt
atcacagcctgtacgtgaaggtgaagggaatgtgttcaaaaacaagcggattctcatggaacacatccacaagctgaaggcagacaaaggcccgcaagaagctcctggctgaccaggctgaggcccgaggtctaagaccaaggaagcacgcaagcggcgtgaagagcgcctcca
ggccaagaaggaggatcatcaagactttatccaagggaagagaccaagaataaa