

IGF2-insulin-like growth factor 2 (somatomedin A) Gene

Reference: PTXBC000531

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C11orf43; GRDF; IGF-II; PP9974; insulin-like growth factor II; T3M-11-derived growth factor; insulin-like growth factor 2 (somatomedin A); insulin-like growth factor type 2; preptin; insulin like growth factor 2		

For research use only.

GENE DESCRIPTION	insulin-like growth factor 2 (somatomedin A)
GENE ID	3481
ORIGIN SPECIES	Human
ACCESSIONS	BC000531
VECTOR	pDONR223
NCBI SYMBOL	IGF2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggaatccaatggggaagtcgatgtgtgtctcactcttctggccttcgctcgtgctgattgctgttaccgccccagtgagaccctgtgcggcggggagctgtgggacaccctccagttcgtctgtggggaccgcggtcttactcagcaggcccgcaagccgtgtgagccg
tcgcagccgtggcatcgttgaggagtgctgttccgagctgtgacctggcctctggagacgtactgtctacccccccaagtcgagagggagctgtgcaccctccgacctgtctccggacaactccccagataccccgtgggcaagtcttccaatagacacctggaagcag
tccaccagcgctgcgaggggctgctgcccctctgctgcccgcgggggtcacgtgctgcgaaggagctcgagggcgttcaggaggccaaacgtaccgtccctgattgcttaccaccaagaccccgccacggggggcggcccccagatggccagcaatcgga
agtgga