

MGLL-monoglyceride lipase Gene

Reference: PTXBC000551

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HU-K5; HUK5; MAGL; MGL; lysophospholipase homolog; monoacylglycerol lipase		
For research use only.			

GENE DESCRIPTION	monoglyceride lipase
GENE ID	11343
ORIGIN SPECIES	Human
ACCESSIONS	BC000551
VECTOR	pDONR223
NCBI SYMBOL	MGLL
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaaacaggacctgaagaccctccagcatgccagaggaaagtccccaggcggaccccgagagcattccctaccaggacctccctcactggcgaatgcagacggacagctctctcaggtactgaaacccacaggcacaccaaggccctcatcttctgtccatgg
agccggagagcagctggccgctatgaagagctggctcgatgctggtggcctggacctgctggtgttcgccacgacctgttggccacggacagagcgaaggaggagagatgtagtctgactccacgttttctcagggatgtgtgcagcatgtggattccatgcagaaag
actaccctggcctctctctctctggtggccactccatgggagcgcacatccctcagggccgcagagagggccgacctcgcggcatggtactcattcgcctctgggtcttccaatctgaatctgcaacaacttcaaggtccttctgctgcaaaagtgtcaacctgtgctg
ccaaactgtccctggcccatcgactccagcgtgctctcctggaataagacagaggtcgacatttataactcagacccctgactcggggcaggcgtgaagggtgtgtctcggcatccaactgctgaatgccgtctcaggggtgagcgcgcctcccaagctgactgtgccctct
gctgctccaggcctctgcccctatgtgacagcaaggggcctacctctcatggagttagccaagagccaggacaagactcgaattatgaaggtgacctaccatgttctccacaaggagcttctgaagtcaccaactccgtctccatgaataaacatgtgggtctctcaag
gacagccacggcaggaaactgctccccacctga