

RABL4-RAB, member of RAS oncogene family-like 4 Gene

Reference: PTXBC000566

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RABL4; BBS19; RAYL; intraflagellar transport protein 27 homolog; RAB, member of RAS oncogene family-like 4; rab-like protein 4; intraflagellar transport 27		

For research use only.

GENE DESCRIPTION	RAB, member of RAS oncogene family-like 4
GENE ID	11020
ORIGIN SPECIES	Human
ACCESSIONS	BC000566
VECTOR	pDONR223
NCBI SYMBOL	RABL4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtgaagctggcagccaaatgcatcctggcaggagaccagcagtggtggaagaccgcccctggcacagatctccgagtgatggagccatttcagaaaaagctacaccctgacacaggaatggatttggtggaagacagtgccagttcctgacacgggagacagtggtgaac
tcttcattttgactctgctgcaaggagctgttttcggaaatgctggataaattgtgggagagtcaccaatgtcttatgtctctctatgatgtgaccaatgaagaatcctcaacaactgcagcaagtgctggagaaggctcggtcacaggctcaggcatctctcccaggtgttttagttggg
aacaagacagacctggccggcagacgagcagtggtgactcagctgagggccgggcatgggctgggcccaggccctggaatgtttgaacatccgtgaaagagatggaagaaactcgaagccctttccactgccttgccaagcagttccaccagctgtaccgggagaagggtgaggttt
tccggccctggcatga