



Reference: PTXBC000574

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PCPE; PCPE-1; PCPE1; procollagen C-endopeptidase enhancer 1; procollagen C-proteinase enhancer 1; procollagen COOH-terminal proteinase enhancer 1; procollagen, type 1, COOH-terminal proteinase enhancer; type 1 procollagen C-proteinase enhancer protein; ty		

For research use only.

GENE DESCRIPTION	procollagen C-endopeptidase enhancer
GENE ID	5118
ORIGIN SPECIES	Human
ACCESSIONS	BC000574
VECTOR	pDONR223
NCBI SYMBOL	PCOLCE
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr