

RPL23AP7-ribosomal protein L23a pseudogene 7 Gene

Reference: PTXBC000596

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	ribosomal protein L23a pseudogene 7
GENE ID	118433
ORIGIN SPECIES	Human
ACCESSIONS	BC000596
VECTOR	pDONR223
NCBI SYMBOL	RPL23AP7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcactcacttcaggcgcccaagacactgcgactccgaggcagccagatcctcagaagacaccccccaggagaacaagcttgccactatgctatcatcaagttccgctgacctgagtcggccgtgaagaagatagaagaaaacaacacgttgttgcactgtggat
gttaagccaacaagcaccagatcagacaggctgtgaagaagctctatgacagtgtgtgccaaggtcaccacccgtgattgtctgataaagagaagaagcatatgttcgactgtctcctgattatgatgctttcaatgtgttaacaaaattggatcacctaaactgagtcagctggcta
a