

DDA1-DET1 and DDB1 associated 1 Gene

Reference: PTXBC000615

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C19orf58; PCIA1; DET1- and DDB1-associated protein 1; PCIA-1; cross-immune reaction antigen PCIA1; placenta cross-immune reaction antigen 1; DET1 and DDB1 associated 1		

For research use only.

GENE DESCRIPTION	DET1 and DDB1 associated 1
GENE ID	79016
ORIGIN SPECIES	Human
ACCESSIONS	BC000615
VECTOR	pDONR223
NCBI SYMBOL	DDA1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcagatttttgaaggactgctgtctacaacaaagcaatttttagtcgatttcacgcggactccgtgtgcaaaagcctgaaccgacggccctcagtcctacccgcgagtcacccgtctgaacagatcatcgtgacagaaaagacaacatcctcctgcgtacctgcacag
caatgggacaaaaaagacgctgccaagaagagagaccaggagcgaatggagctggaaggcgagagctccgacctcccccaaggtggcgcgacccagacagccagacatgcacgaggacacttaa